




New ESS Co., Ltd.



Order NO. F070S
ISO NO.500.104 274.175.070

色码
color coding

出厂编号 (Order NO.): F070S
Burs snape(A)
磨头形状

磨头的齿型是中等单纹
Standard the burs tooth number

磨头直径(D)7mm
Burs diameter(D)7mm

国际编号: (ISO NO.) .500 104 274 175 070

磨头材质、硬质合金
Burs material, Carbide


直径 (D)7mm
Diameter(D)7mm

齿数: 标准单刃
Tooth number: standard single blade


柄的直径/长度 (2.35mm/44.5mm)
Shank diameter/length(2.35mm/44.5mm)

磨头形状 (A) Burs shape(A)

Cut, medium, (ISO-编号, 175), bule ring, cut carbide effectively and create smooth effect.



Order NO.	A0613S	C016S	L023S	D023S	D030S	D040S	D050S	C060S	E060S	F060S	F070S	M060S
ISO NO.500.104	107.175.060	141.175.016	194.175.023	001.175.023	001.175.030	001.175.040	001.175.050	144.175.060	277.175.060	274.175.060	274.175.070	161.175.060



Order NO.	LF060S	LF070S	CN050S	CN060S	G060S	L040S	L060S	L060S	L070S	H040S	MA023S
ISO NO.500.104	263.175.060	263.175.070	239.175.050	239.175.060	257.175.060	194.175.040	194.175.045	194.175.060	194.175.070	275.175.040	187.175.023

Tel: +86-400-053-0087

Fax: +86-400-676-1100/1247


Email: info@dentalburs.cn

Web: http://www.dentalburs.cn




New ESS Co., Ltd.


Crosscut, medium, (ISO-NO.190). blue ring, apply to cut extensive alloy and acrylics.



Cut, coarse, (ISO-NO.215), green ring, cut acrylics, plaster and such nonmetal material rapidly and effectively.



Crosscut, coarse, (ISO-NO.220), green ring, apply to cut extensive alloy and acrylics denture and plaster.



Tel: +86-400-053-0087

Email: info@dentalburs.cn




Fax: +86-400-676-1100/1247

Web: <http://www.dentalburs.cn>




New ESS Co., Ltd.


Crosscut "X", fine, (ISO-NO.140), red ring, apply to the fragile and hard-cutting carbide, acrylics and porcelain.



Crosscut "X", superfine, (ISO-NO.110), yellow ring, apply to the final finishing of porcelain, noble metal and alloy.



Carbide bur covered by titanium layer



Tel: +86-400-053-0087

Email: info@dentalburs.cn

Fax: +86-400-676-1100/1247

Web: http://www.dentalburs.cn