




# NEW ESS CO.,LTD.


| Item        | Shank Diameter | Length | ●   | ●  |
|-------------|----------------|--------|---|--|
| FG          | 1.6mm          | 19mm   | Material Size:diameter 1.6mm/length 19mm    | Material Size: diameter 2.35mm/length 19mm |
| FGOS        | 1.6mm          | 21mm   | Material Size:diameter 1.6mm/length 21mm    | Material Size: diameter 2.35mm/length 26mm |
| FG Surgical | 1.6mm          | 25mm   | Material Size:diameter 1.6mm/length 25mm    | Material Size: diameter 2.35mm/length 44mm |
| RA          | 2.35mm         | 22.5mm | Material Size:diameter 2.35mm/length 22.5mm | /  |
| RA Surgical | 2.35mm         | 26mm   | Material Size:diameter 2.35mm/length 26mm   | /  |
| HP          | 2.35mm         | 44mm   | Material Size:diameter 2.35mm/length 44mm   | /  |


|  |             |         |         |         |         |         |         |         |         |         |         |   |
|--|-------------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---|
|  平头柱状裂钻<br>Cross-Cut Cylinder | Order No.   | 556     | 557     | 558     | 559     | 560     | 561     | 562     | 563     | 557L    | 558L    |   |
|  | ISO NO.     | 107/009 | 107/010 | 107/012 | 107/014 | 107/016 | 107/018 | 107/021 | 107/023 | 107/010 | 107/102 |   |
|  | FG          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | /       | / |
|  | FGOS        | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ● |
|  | FG Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ● |
|  | RA          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ● |
|  | RA Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ● |
|  | HP          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ● |


|  |             |         |         |         |         |         |         |         |         |         |         |
|--|-------------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|
|  球钻<br>Round | Order No.   | 1/4     | 1/2     | 1       | 2       | 3       | 4       | 5       | 6       | 7       | 8       |
|  | ISO NO.     | 001/005 | 001/006 | 001/008 | 001/010 | 001/012 | 001/014 | 001/016 | 001/018 | 001/021 | 001/023 |
|  | FG          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | FGOS        | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | FG Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | RA          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | RA Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | HP          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |

|   |             |         |         |         |         |         |         |         |         |         |         |
|---|-------------|---------|---------|---------|---------|---------|---------|---------|---------|---------|---------|
|  平头裂钻<br>Cross-cut Taper | Order No.   | 698     | 699     | 700     | 701     | 702     | 703     | 704     | 700L    | 701L    | 702L    |
|   | ISO NO.     | 168/008 | 168/009 | 168/010 | 168/012 | 168/016 | 168/021 | 168/023 | 171/010 | 171/012 | 171/016 |
|   | FG          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | /       | /       | /       |
|   | FGOS        | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|   | FG Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|   | RA          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|   | RA Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|   | HP          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |


# NEW ESS CO.,LTD.


|  |             |         |         |         |         |         |         |         |         |
|--|-------------|---------|---------|---------|---------|---------|---------|---------|---------|
|  平头柱状<br>Cylinder | Order No.   | 56      | 57      | 58      | 59      | 60      | 57L     | 58L     | 60L     |
|  | ISO NO.     | 107/009 | 107/010 | 107/012 | 107/014 | 107/016 | 107/010 | 107/012 | 107/016 |
|  | FG          | ●       | ●       | ●       | ●       | ●       | /       | /       | /       |
|  | FGOS        | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | FG Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | RA          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | RA Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
| HP   | ●           | ●       | ●       | ●       | ●       | ●       | ●       | ●       |         |


|   |             |         |         |         |         |         |         |         |         |
|---|-------------|---------|---------|---------|---------|---------|---------|---------|---------|
|  倒锥<br>Inverted cone | Order No.   | 34      | 35      | 36      | 37      | 38      | 39      | 40      | 41      |
|   | ISO NO.     | 010/008 | 010/010 | 010/012 | 010/014 | 010/016 | 010/018 | 010/021 | 010/023 |
|   | FG          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|   | FGOS        | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|   | FG Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|   | RA          | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|   | RA Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
| HP  | ●           | ●       | ●       | ●       | ●       | ●       | ●       | ●       |         |


|  |             |         |         |         |         |         |         |         |
|--|-------------|---------|---------|---------|---------|---------|---------|---------|
|  平头锥形<br>Taper | Order No.   | 168     | 169     | 170     | 171     | 172     | 170L    | 171L    |
|  | ISO NO.     | 168/008 | 168/009 | 168/010 | 168/012 | 168/016 | 171/010 | 171/012 |
|  | FG          | ●       | ●       | ●       | ●       | ●       | /       | /       |
|  | FGOS        | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | FG Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | RA          | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | RA Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
| HP   | ●           | ●       | ●       | ●       | ●       | ●       | ●       |         |


# NEW ESS CO.,LTD.

|  |             |         |         |         |         |         |         |         |
|--|-------------|---------|---------|---------|---------|---------|---------|---------|
|  <p>梨形<br/>Pear</p> | Order No.   | 329     | 330     | 331     | 332     | 331L    | 332L    | 333L    |
|  | ISO NO.     | 232/006 | 232/008 | 232/010 | 232/012 | 232/010 | 232/012 | 232/014 |
|  | FG          | ●       | ●       | ●       | ●       | ●       | /       | /       |
|  | FGOS        | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | FG Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | RA          | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
|  | RA Surgical | ●       | ●       | ●       | ●       | ●       | ●       | ●       |
| HP   | ●           | ●       | ●       | ●       | ●       | ●       | ●       |         |

|   |             |         |         |         |         |         |         |
|---|-------------|---------|---------|---------|---------|---------|---------|
|  <p>圆头锥状<br/>Round End Taper</p> | Order No.   | 1169    | 1170    | 1171    | 1172    | 1171L   | 1172L   |
|   | ISO NO.     | 199/009 | 199/010 | 199/012 | 199/016 | 199/012 | 199/016 |
|   | FG          | ●       | ●       | ●       | ●       | /       | /       |
|   | FGOS        | ●       | ●       | ●       | ●       | ●       | ●       |
|   | FG Surgical | ●       | ●       | ●       | ●       | ●       | ●       |
|   | RA          | ●       | ●       | ●       | ●       | ●       | ●       |
|   | RA Surgical | ●       | ●       | ●       | ●       | ●       | ●       |
| HP  | ●           | ●       | ●       | ●       | ●       | ●       |         |

|   |             |         |         |         |         |         |
|---|-------------|---------|---------|---------|---------|---------|
|  <p>圆头柱状<br/>Cylinder</p> | Order No.   | 1155    | 1156    | 1157    | 1158    | 1159    |
|   | ISO NO.     | 137/008 | 137/009 | 137/010 | 137/012 | 137/014 |
|   | FG          | ●       | ●       | ●       | ●       | ●       |
|   | FGOS        | ●       | ●       | ●       | ●       | ●       |
|   | FG Surgical | ●       | ●       | ●       | ●       | ●       |
|   | RA          | ●       | ●       | ●       | ●       | ●       |
|   | RA Surgical | ●       | ●       | ●       | ●       | ●       |
| HP  | ●           | ●       | ●       | ●       | ●       |         |

|   |             |         |         |         |         |
|---|-------------|---------|---------|---------|---------|
|  <p>圆头裂钻<br/>Round End Cross-cut Taper</p> | Order No.   | 1701    | 1702    | 1703    | 1704    |
|   | ISO NO.     | 199/012 | 199/016 | 199/021 | 199/023 |
|   | FG          | ●       | ●       | ●       | ●       |
|   | FGOS        | ●       | ●       | ●       | ●       |
|   | FG Surgical | ●       | ●       | ●       | ●       |
|   | RA          | ●       | ●       | ●       | ●       |
|   | RA Surgical | ●       | ●       | ●       | ●       |
| HP  | ●           | ●       | ●       | ●       |         |

|  |             |         |         |         |
|--|-------------|---------|---------|---------|
|  <p>圆头柱状列钻<br/>Round End Corss-cut Cylinder</p> | Order No.   | 1556    | 1557    | 1557    |
|  | ISO NO.     | 137/009 | 137/010 | 137/012 |
|  | FG          | ●       | ●       | ●       |
|  | FGOS        | ●       | ●       | ●       |
|  | FG Surgical | ●       | ●       | ●       |
|  | RA          | ●       | ●       | ●       |
|  | RA Surgical | ●       | ●       | ●       |
| HP   | ●           | ●       | ●       |         |